

Post-Op Knee Bracing

DonJoy Telescoping TROM



Ideal Application

TROM offers state-of-the-art range-of-motion control and locking capabilities. Only TROM allows flexion and extension adjustments at the touch of a button and a quick lock from 0° to 90° in 10° increments.

Features and Benefits

- **NEW! Telescoping bar**
- Innovative, easy-to-use hinge
- Quick-lock mechanism from 0° to 90°
- Easily size and custom fit the brace to the patient
- Buckles allow for fast and easy strapping

Part Number:

11-0961-9-06000	Enhanced TROM standard
11-0962-9-06000	Enhanced TROM, Shells
11-0964-9-06000	Enhanced TROM, XL

DonJoy TROM Advanced Cool



Ideal Application

TROM Advance offers range of motion control and locking capabilities following ACL, PCL, MCL, and LCL surgeries, Meniscal repairs, Patella realignment, Regenerative Chondroplasty, Stable femoral fractures, Total Knee replacements and High Tibial Osteotomies

Features and Benefits

- Tele-Fit: Sliding mechanism allowing independent strap movement just proximal and distal of the knee to accommodate cold therapy wraps and dressings and for increased strap adjustability

Part Number:

11-9114-9	TROM Advance, Cool
11-9116-9	TROM Advance, XL Cool

DonJoy Telescoping IROM



Ideal Application

The IROM is the industry standard providing precise range-of-motion control in 10 degree increments. The IROM is ideal for ACL, PCL, MCL, and LCL surgeries. The telescoping bars allow the brace to fit a variety of patients.

Features and Benefits

- **NEW! Telescoping bar**
- Precise range of motion control
- Easily size and custom fit the brace to the patient
- Buckles allow for fast and easy strapping
- Brace can extend to ankle to help prevent migration

Part Number:

11-0968-9-13066	Enhanced IROM, Standard
11-0971-9-13066	Enhanced IROM, XL

Accessories:

11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left

DonJoy LS Bilateral w/ & w/o without Shells



Ideal Application

The short length of the ELS Bilateral facilitates rehab exercises and can provide full extension support for rheumatoid arthritis or diabetic neuropathy patients. The ELS Bilateral is available with or without shells.

Features and Benefits

- Lock at full extension
- Buckles allow for fast and easy strapping
- 50° flexion stop kit available for men

Part Number:

11-2128-9-13066	Enhanced ELS, Bilateral
11-2129-9-13066	Enhanced, ELS, Bilateral, Shells

DonJoy Telescoping Cool ELS



Ideal Application

The ELS features a slide lock mechanism for the ability to lock at full extension or have free range of motion. The brace is very easy to use. The brace is ideal for post-op use when full extension weight-bearing protocols are used. The Cool ELS gives improved patient comfort with abbreviated foam.

Features and Benefits

- **NEW! Telescoping bar**
- Lock at full extension
- Easily size and custom fit the brace to the patient
- Brace can extend to ankle to help prevent migration


Part Number: 11-0967-9-13066 Enhanced ELS, Cool

Accessories

11-0080-9-00000	50° flexion stop, ELS
11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left

Post-Op Knee Bracing

DonJoy Telescoping



Ideal Application

The ELS features a slide lock mechanism providing the ability to lock at full extension or have free range of motion. The brace is very easy to use with a simple yet effective design. The brace is ideal for post-op use when full extension weight-bearing protocols are used.

Features and Benefits

- **NEW!** Telescoping bar
- Lock at full extension
- Easily size and custom fit the brace to the patient
- Brace can extend to ankle to help prevent migration
- 50° flexion stop kit available for meniscal repair protocols.

Part Number: 11-0965-9-13066 Enhanced ELS, Standard

Accessories:

11-0080-9-00000	50 flexion stop, ELS
11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left

DonJoy Cool IROM



Ideal Application

For warmer climates, the Cool IROM provides immobilization and precise range-of-motion control for ACL, PCL, MCL and LCL surgeries. Formable hinge bars allow for varus/valgus positioning

Features and Benefits

- 10° incremental settings
- 4 straps and abbreviated foam at key control points
- Integrated strap/hinge bar design
- Hinge bar breakdown capability

Part Number:

11-0175-2-13066	Cool IROM	Short (<32"/80cm)
11-0175-4-13066	Cool IROM	Long (>32"/80cm)

Accessories:

11-0024-x	Ankle stirrup, right
11-0026-x	Ankle stirrup, left

DonJoy PCL Fource®



Ideal Application

For added stiffness and tissue containment, the IROM with shells offers full-circumference shells. The knee brace offers effective post-operative immobilization and protected range of motion associated with ligament and fracture related surgeries

Features and Benefits

- Reverse Fourcepoint hinge technology
- Tele-Fit allows independent strap movement proximal and distal of knee allowing an overall better fit
- Telescoping adjustability from 17 3/4" to 24"
- Quick release buckles

Part Number:

11-0488	PCL Fource	< 26"
11-0489	PCL Fource, XL	26" - 32"

Measure 6" above mid-patella

DonJoy Rehab TROM



Ideal Application

Ideal for total knees, meniscus repairs, regenerative chondroplasty, ligament surgeries and patella realignments, the Rehab TROM gives range of motion control and short cuffs for extended rehabilitation. It can be used as a rehab brace or the cuffs can be purchased separately to be used with the TROM post-op brace.

Features and Benefits

- Innovative, easy-to-use hinge
- Quick release buckles
- Removable, washable, latex-free cuffs
- Easy to apply

Part Number:

11-0295-x-06000	Rehab TROM
11-0745-x-06000	TROM Conversion Kit

Size Information

x=2=S	14"-18" (36-46cm)	x=4=L	22"-27" (56-69cm)
x=3=M	18"-22" (46-56cm)	x=5=XL	27"-32" (69-81cm)

DonJoy TROM Adjuster



Ideal Application

Create the optimal knee off-load with just a few turns. TROM Adjuster is the first post-op brace to offer controlled and consistent off-loading to create a better healing environment. TROM Adjuster represents a comprehensive post operative solution for meniscal and chondral defect procedures

Features and Benefits


- Off-loading adjustable hinge technology
- Tele-Fit: Sliding mechanism allowing independent strap movement just proximal
- Telescoping adjustability between 17 3/4" and 24"
- Accommodates cold therapy wraps/dressings
- Sili grip strap pads assist in increasing brace

Part Number:

11-5117-9	TROM Adjuster XL, Medial RT/Lateral LT	26-32"
11-5118-9	TROM Adjuster XL, Medial LT/Lateral RT	26-32"
11-5115-9	TROM Adjuster, Medial R/Lateral LT	< 26"
11-5116-9	TROM Adjuster, Medial LT/Lateral RT	< 26"

Post-Op Knee Bracing

IROM Bilateral with and without Shells



Ideal Application

The IROM is the industry standard providing precise range-of-motion control in 10 degree increments. The IROM is ideal for ACL, PCL, MCL, and LCL surgeries. The telescoping bars allow the brace to fit a variety of patients.

Features and Benefits

- 10 degree incremental range of motion control
- Quick-release buckles allow for fast and easy strapping
- Shells provide added stability

Part Number:

11-0972-9-13066 Enhanced IROM Bilateral w/ shells
11-2127-9-13066 Enhanced IROM Shells

Breg T-SCOPE



Ideal Application

The T-Scope features telescoping hinge bars that extend from 17" - 27", along with a convenient drop-lock style hinge that controls range of motion 10° increments. The T-Scope also features friction locks to easily secure the brace in position once the desired length is selected, along with quick clip buckles and strap lock clips for easy application and fitting.

Features and Benefits

Telescoping hinge bars extendable from 17" – 27" (43.2 cm - 68.6 cm) accommodates a larger range of patient leg circumference). XL version available *(high circumference 27" – 35"). Lightweight - Weighs less than 2.5 lbs. Strap lock clips keep straps positioned on bars. Quick clip buckle easy on/off. Quick lock offers ability to lock out brace from –10 to 30°. Reduced foam for greater breathability. Full foam version available.

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Part #	Description
07714	T-Scope Post-Op
07715	T-Scope Post-Op, XL
07716	T-Scope Post-Op, Full Foam
07717	T-Scope Post-Op, Full Foam, XL

Breg Post-Op



Featuring the premier T-Scope hinge, the upper and lower hinge bars on the Post-Op can be easily contoured or shortened for an optimized fit, and the foam can be trimmed to suit smaller patients.

Product Features and Benefits

Reduced foam is perfect for warmer climate. Spring loaded flexion/extension stops for easy range of motion control from 0° to 120°. Ability to lock-out in full extension . Scored anodized bars for easy breakdown to rehab style brace. Available in short and long. Available in XL versions. Latex-free. Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Size and Part

Part #	Description
00112	Post-Op, Short
00113	Post-Op, Short, XL
00114	Post-Op, Long
00115	Post-Op, Long, XL

Note: Short <32" inseam. Long >32" inseam. XL = 27" - 31" thigh circumference (measurement taken 6" above mid-patella).

Breg Post-Op Lite



Featuring the premier T-Scope hinge, the Post-Op Lite incorporates a drop-lock style hinge that allows the user to easily control the range of motion in 10° increments. The upper and lower hinge bars on the Post-Op can be easily contoured or shortened for an optimized fit, and it has reduced foam coverage for greater breathability and comfort.

Product Features and Benefits

Telescoping hinge bars extendable from 17" – 27" (43.2 cm - 68.6 cm) accommodates a larger range of patient leg lengths. Universally sized (fits up to 30.5" thigh circumference). XL version available *(high circumference 27" – 35"). Lightweight - Weighs less than 2.5 lbs. Strap lock clips keep straps positioned on bars. Quick clip buckle easy on/off. Quick lock offers ability to lock out brace from –10 to 30°. Reduced foam for greater breathability. Full foam version available.

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Size and Part

Part #	Description
00116	Post-Op Lite, Short
00117	Post-Op Lite, Short, XL
00118	Post-Op Lite, Long
00119	Post-Op Lite, Long, XL

Note: Short <32" inseam. Long >32" inseam. XL = 27" - 31" thigh circumference (measurement taken 6" above mid-patella).

Post-Op Knee Bracing

Breg Post-Op Rehab



The Post-Op Rehab brace is a shortened version of BREG's Post-Op brace, measuring 15" in overall length. This version is commonly used for more progressive protocols while assisting patients to return to normal activities.

Product Features and Benefits

Anodized hinge bars. Spring loaded flexion/extension stops for easy range of motion control from 0° to 120°. Full foam for comfort and tissue containment. Latex-free

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
0123	Post-Op Rehab

Note: Short Fits up to 27" thigh circumference (measurement taken 6" above mid-patella).

Breg Quick Fit Post-Op



The Quick-Fit Post-Op is pre-sized for ease of application. It combines user friendly Velcro® straps, and is one of the easiest Post-Op braces to apply.

Product Features and Benefits

Spring loaded flexion/extension stops allow for easy range of motion from 0° to 120°. Anodized hinge bars (may not be broken down). Full 1/2" foam provides comfort and soft tissue containment. Pre-sized S through XL for easier application. Quick fit buckles for easy on/off. Four straps for quicker application. Available in short and long. Latex-free.

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
0014 X	Quick Fit Post-Op, Short, S-XL
0015X	Quick Fit Post-Op, Long, S-XL

Size Information

SIZE #	SIZE	THIGH CIRCUMFERENCE
X=2	S	14-18" (36-46 CM)
X=3	M	18-22" (46-56 CM)
X=4	L	22-27" (56-69 CM)
X=5	XL	27-31" (69-79CM)

Breg EPO Lite



An economical solution for post-operative indications, the EPO Lite features the Quick-Set hinge and reduced foam coverage for greater breathability and comfort.

Product Features and Benefits

Easy-to-move set pins for range of motion. Ability to lock-out in full extension. Scored bars for easy breakdown to rehab style brace. Available in short and long. Available in a lite version. Available in a quick fit version.

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
09360	Economy Post-Op Lite, Short 23.5" lgth
09370	Economy Post-Op Lite, Long 26.5" lgth

Note: Short <32" inseam. Long >32"

Post-Op Knee Bracing

Breg Post-Op with Shells



For additional rigidity and stability, the shells on the Post-Op with Shells enable full circumferential tissue containment. The polyurethane shells are adjustable and universally sized.

Product Features and Benefits

Full circumferential shells for tissue containment and extra support. Spring-loaded flexion/extension stops allow for easy range of motion from 0° to 120°. Ability to lock-out in full extension. Available in short and long. Latex-free

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
00121	Post-Op w/ Shells, Short (23.5" length)
00122	Post-Op w/ Shells, Long (26.5" length)

Breg EPO



The EPO is an economical alternative to our premier Post-Op line. EPOs feature a Quick-Set hinge providing complete range of motion control.

Product Features and Benefits

Easily adjustable set pins for range of motion control. Ability to lock-out in full extension, as well as a range of flexion/extension settings. Scored bars for easy breakdown to rehab style brace. Available in short and long. Latex-free

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
09340	Economy Post-Op, Short (23.5" length)
09350	Economy Post-Op, Long (26.5" length)

Note: Short <32" inseam. Long >32"

Bledsoe Lever Lock Knee Brace



The **Bledsoe Lever Lock Knee Brace** is designed to provide a quick extension lock of the arthroscopic ACL reconstruction. The levers quickly unlock to allow motion for rehabilitation or passive motion exercise, but may be quickly locked "rigidly" in full extension for ambulation or sleeping. Application of this brace is the easiest on the market.

Product Features and Benefits

Hinge upright pockets attach quickly to foam wraps with hook-to-pile closures. One piece strap-to-wrap for faster attachment. A small anti-migration strap is laced through the foam wrap to reduce slippage of the brace. The "regular" or "long" sizes fit left or right legs of most patients. Extra-long version available for inseams of 32" (81 cm) or greater. Tamper-resistant tie seals for hinge are provided for patient compliance.

Indications: immediate application following arthroscopic reconstruction of the anterior cruciate ligament or other procedures or injuries which require locking the knee in full extension with intermediate periods of full range of motion.

Size Information and Order Information

Brace Size	Part Number	Replacement Cuffs
Regular	AK070005	AK073005
Long	AK070007	AK073007
X-Long	AK070009	AK073009

Post-Op Knee Bracing

Breg Quick Fit EPO



The Quick Fit EPO is pre-sized for ease of application in the OR and recovery room.

Product Features and Benefits

Easy-to-move set pins for range of motion. Ability to lock-out in full extension. Full foam wrap for comfort and compression. Easy to use Velcro® straps. Available in short and long.

Indications: ACL, PCL, MCL, LCL injuries. Tibial plateau fractures. Osteochondral repairs. Meniscal repairs. Patella tendon repairs. Condylar fractures. Sprains/strains of the knee. Total knees. HTOs (high tibial osteotomy).

Order Information

Part #	Description
0934X	Quick Fit EPO, Short (23.5" length), S-XL
0935X	Quick Fit EPO, Long (26.5" length), S-XL

Size Information

SIZE #	SIZE	THIGH CIRCUMFERENCE
X=1	XS	12-15" (30-38 CM)
X=2	S	15-18" (38-46 CM)
X=3	M	18-21" (46-53 CM)
X=4	L	21-24" (53-61 CM)
X=5	XL	24-27" (61-69 CM)
X=6	XXL	27-30" (69-79 CM)

Note: Circumference taken 6" above mid patella

Bledsoe Merit OR Knee Brace



The **Bledsoe Merit OR Knee Brace** frees the surgeon or authorized technician from the time-consuming post-op fitting and application that is required with more complex brace designs. The **Bledsoe Cool Merit OR** knee brace provides the same benefits as the Bledsoe OR in a cool, open design for warm climate.

Product Features and Benefits

Fits any patient in just four sizes. Vinyl hinge upright pockets with hook-to-pile closures for quick attachment to wraps. Tamper-resistant tie seals supplied for patient security and compliance. NEW! **Post Op Hinge** provides fast and easy adjustment with flexion limits of -10 to 110° and extension limits of -10 to 110° in 10° increments. No tools to adjust. The hinge can be locked in position from -10 to 40°. Can be removed in one piece for easy reattachment.

Indications: This device is indicated for locked or limited motion control of knee during rehabilitation after operative procedures or injury to knee ligaments, cartilage, or stable or internally fixed fractures of tibial plateau, condyles, or proximal tibia and distal femur.

Size and Order Information

Length	Brace Size	Merit OR	Cool Merit OR
16-22" (40.6-55.9 cm)	Short	EK080001	EK083001
22-26" (55.9-66 cm)	Regular	EK080005	EK083005
26-30" (66-76.2 cm)	Long	EK080007	EK083007
30"+ (76.2 cm+)	X-Long	EK080009	EK083009

Post-Op Knee Bracing

Bledsoe Thruster RLF



The **Bledsoe Thruster RLF OA Knee Brace** uses muscle power to place a medially (Thruster RLF Medial) or laterally (Thruster RLF Lateral) directed force against the knee during terminal extension to help relieve pain for moderate to severe symptomatic lateral or medial compartment arthrosis.

Indications: *Thruster RLF LA* - This device is indicated for Lateral Unicompartmental Gonarthrosis, to help relieve pain and inflammation, to aid in delaying osteotomy or total joint procedures, to increase daily walking time while increasing activities, and to aid in slowing the progression of the condition. Also indicated for postoperative or post-injury applications following medial collateral ligament injury or reconstruction to help relieve stress on these healing structures.

Thruster RLF MA - This device is indicated for Medial Unicompartmental Gonarthrosis, to help relieve pain and inflammation, to aid in delaying osteotomy or total joint procedures, to increase daily walking time while increasing activities, and to aid in slowing the progression of the condition. Also indicated for postoperative or post-injury applications following lateral and postero-lateral injuries to the connective tissues of the knee to help relieve stress on these healing structures.

Size Information and Order Information

Thigh Circumference	Calf Circumference		Men				Women	
			Standard		Short		Standard	
Men and Women	Men	Women	MA	LA	MA	LA	MA	LA
15"-17.5" (38cm-44.5cm)	12"-13.5" (30.5cm-34cm)	10.5"-12" (26.7cm-30.5cm)	XS	XS	XS	XS	XS	XS
17.5"-20" (44.5cm-51cm)	13.5"-15" (34.5cm-38cm)	12"-13.5" (30.5cm-34cm)	S	S	S	S	S	S
20"-22.5" (51cm-57cm)	15"-16.5" (38cm-41.9cm)	13.5"-15" (34cm-38cm)	M	M	M	M	M	M
22.5"-25" (57cm-63.5cm)	16.5"-18" (41.9cm-45.7cm)	15"-16.5" (38cm-41.9cm)	L	L	L	L	L	L
25"-27.5" (63.5cm-70cm)	18"-19.5" (45.7cm-49.5cm)	16.5"-18" (41.9cm-45.7cm)	XL	XL	XL	XL	XL	XL

1. Measure leg circumference 6" (15 cm) above and below the femoral epicondyle.
2. Note patient height, sex, and brace type in table.

Please contact Customer Service for order information

Bledsoe Crossover ROM



The **Bledsoe Crossover ROM** post-operative knee brace is designed for very mild meniscal and/or collateral injuries. Patients with minor knee ligament strains and stretches or meniscal damage can benefit from added support to prevent more serious injury. The Crossover ROM provides control over range-of-motion of the knee to aid in healing and rehabilitation.

Product Features and Benefits

Cool, breathable Spacer Mesh provides support, compression and warmth to the knee. Easily adjustable ROM hinges for ideal anatomical placement. Front closure design for easy application. Interchangeable Sport Hinges convert the Crossover ROM to the Crossover.

Indications: locked or limited motion control during rehabilitation after arthroscopic procedures or injury.

Size and Order Information

Thigh Circumference	Brace Size	Crossover ROM	Sport Hinges*
15"-17.5" (38.1-44.5 cm)	XS	RK090001	RK066301
17.5"-20" (44.5-50.8 cm)	S	RK090003	RK066303
20"-22.5" (50.8-57.2 cm)	M	RK090005	RK066305
22.5"-25" (57.2-63.5 cm)	L	RK090007	RK066307
25"-27.5" (63.5-69.9 cm)	XL	RK090009	RK066309
27.5"-29" (69.9-73.7 cm)	2XL	RK090011	RK066311

Post-Op Knee Bracing

Bledsoe Original Knee Brace



The **Bledsoe Original Knee Brace** gives fracture treatment patients maximum support and comfort. The Original, our most adaptable brace, integrates fracture plates as well as shoe inserts, boots and/or pelvic brace for KAFO and/or HKAFO.

Product Features and Benefits

Greatest control and adaptability. Applies in 5 minutes. The Bledsoe Original Knee Brace saves 20 to 30 minutes over traditional cast bracing methods. Offering cast-like support with a weight under 4 lbs. (1.8 kg), the brace gives support in all planes while allowing muscles to maintain their natural shape. Rigid medial and lateral stays and molded posterior supports provide a thorough support system. Non-stretch straps and breathable foam wraps help ensure patient comfort. NEW! **Post Op Hinge** provides fast and easy adjustment with flexion limits of -10 to 110° and extension limits of -10 to 110° in 10° increments. No tools to adjust. The hinge can be locked in position from -10 to 40°. 9 sizes promote best fit. Optional shoe insert provides additional rotational control of the foot and lower leg. Tibial fracture plate combined with the shoe insert convert the Bledsoe Original into a KAFO.

Indications: This device is indicated for immediate application to knees for locked or limited motion control of knee during rehabilitation, after operative procedures or injury to knee ligaments, cartilage, or stable or internally fixed fractures of the tibial plateau, condyles, or proximal tibia and distal femur. Also indicated for immediate post-injury application for stable, non-displaced or internally fixed fractures of the middle to distal femur and the proximal to middle tibia or fibula when using the optional fracture

Size Information and Order Information

Leg Length	Brace Size	Bledsoe Original	Hinge & Thigh Attachment *
17-19" (43-48 cm)	18" (45.7 cm)	EK020018	EK025018
19-21" (48-53 cm)	20" (50.8 cm)	EK020020	EK025020
21-23" (53-57 cm)	22" (55.8 cm)	EK020022	EK025022
23-25" (57-63 cm)	24" (61 cm)	EK020024	EK025024
25-27" (63-69 cm)	26" (66 cm)	EK020026	EK025026
27-29" (69-74 cm)	28" (71 cm)	EK020028	EK025028
29-31" (74-79 cm)	30" (76.2 cm)	EK020030	EK025030
31-33" (79-84 cm)	32" (81.3 cm)	EK020032	EK025032
33-35" (84-89 cm)	34" (86.4 cm)	EK020034	EK025034
35+ (89+ cm)	Contact us for ordering instructions.		

Optional Components

Sz	Shoe Insert		Femoral Plate		Tibial Plate		Winged Tibial Plate	
	Left	Right	Left	Right	Left	Right	Left	Right
S	AM005101	AM005201	AJ005101	AJ005201	AL005101	AL005201	AL004101	AL004201
M	AM005105	AM005205	AJ005105	AJ005205	AL005105	AL005205	AL004105	AL004205
L	AM005109	AM005209	AJ005109	AJ005209	AL005109	AL005209	AL004109	AL004209

Bledsoe Three Panel Knee Immobilizer



The **Three Panel Knee Immobilizer** features: 3/8-inch pile foam laminate with a comfy terry interior; Custom fit up to a 30-inch circumference thigh with the two detachable circumferential panels; Fit and shape to legs with or without dressings using posterior and medial/lateral stays; Patella strap for ultimate knee control.

Indications: for motion control of the knee after malalignment injury, post operative care or during rehabilitation of knee ligaments, cartilage or stable or internally fixed fractures of the tibial plateau, condyles or proximal tibia and distal femur.

Order Information

Length	Product Number
12" (30cm)	SA608012
16" (41cm)	SA608016
20" (51cm)	SA608020
24" (61cm)	SA608024

Post-Op Knee Bracing

Bledsoe Aligner RLF



The **Bledsoe Aligner ESE OA Knee Brace** uses muscle power to place a laterally (Aligner ESE LA) or medially (Aligner ESE MA) directed force against the knee during terminal extension to help relieve pain for mild to moderate symptomatic lateral or medial compartment arthrosis.

Indications: Aligner ESE LA - This device is indicated for Lateral Unicompartmental Gonarthrosis, to help relieve pain and inflammation, to aid in delaying osteotomy or total joint procedures, to increase daily walking time while increasing activities, and to aid in slowing the progression of the condition. Also indicated for postoperative or post-injury applications following medial collateral ligament injury or reconstruction to help relieve stress on these healing structures.

Aligner ESE MA - This device is indicated for Medial Unicompartmental Gonarthrosis, to help relieve pain and inflammation, to aid in delaying osteotomy or total joint procedures, to increase daily walking time while increasing activities, and to aid in slowing the progression of the condition. Also indicated for postoperative or post-injury applications following lateral and postero-lateral injuries to the connective tissues of the knee to help relieve stress on these healing structures.

Size Information and Order Information

Leg Circumference			Brace Size	Aligner MA (Medial)		Aligner LA (Lateral)	
Thigh	Calf (Men)	Calf (Women)		Left	Right	Left	Right
15"-17.5" (38-44.5 cm)	12"-13.5" (30.5-34.5 cm)	10.5"-12" (26.7-30.5cm)	XS	AR079101	AR079201	AR078101	AR078201
17.5"-20" (44.5-51 cm)	13.5"-15" (34.5-38 cm)	12"-13.5" (30.5-34.5 cm)	S	AR079103	AR079203	AR078103	AR078203
20"-22.5" (51-57 cm)	15"-16.5" (38-41.9 cm)	13.5"-15" (34.5-38 cm)	M	AR079105	AR079205	AR078105	AR078205
22.5"-25" (57-63.5 cm)	16.5"-18" (41.9-45.7 cm)	15"-16.5" (38-41.9 cm)	L	AR079107	AR079207	AR078107	AR078207
25"-27.5" (63.5-70 cm)	18"-19.5" (45.7-49.5 cm)	16.5"-18" (41.9-45.7 cm)	XL	AR079109	AR079209	AR078109	AR078209

Universal AFO Attachment: TK074075

Measure the thigh and calf circumferences 6" (15 cm) above and below the adductor tubercle (femoral epicondyle). **NOTE: If either measurement falls outside of the specified ranges, semi-custom versions are available. Specify left or right when ordering. HCFA Code: L1843**

Bledsoe Extender Plus



The **Bledsoe Extender Plus Brace** is Adjustable from 22" to 34" (56 cm to 86 cm) in length and up to 34" (86 cm) in leg circumference, this brace easily fits patients from 5' (152 cm) to 6' 5" (196 cm) tall, and adapts as the patient progresses. The Bledsoe Extender Plus is so versatile there is no need for an expensive custom-fitted brace.

Product Features and Benefits

Fast, easy push-button length adjustment — no tools required. Quickly applies with fast, wrap-around strap design and no D-rings. Cool, open design reduces skin breakdown over long-term use. Four-point attachment and long-hinge arms provide superior control for recovering patients. Stable suspension from formed, double padded ankle plates. **NEW!** Over-center buckles for easy application. **NEW! Post-Op Hinge** provides fast and easy adjustment with flexion limits of -10 to 110° and extension limits of -10 to 110° in 10° increments. No tools to adjust. The hinge can be locked in position from -10 to 40°. Hinge tie seals supplied for patient security and compliance. Optional shoe insert provides additional rotational control of the foot and lower leg.

Indications: This device is indicated for immediate application to knees for locked or limited motion control of knee during rehabilitation after operative procedures or injury to knee ligaments (MCL, LCL, ACL and PCL), patellar ligament, articular cartilage, meniscus, or stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur.

Size Information and Order Information

Size	Part #	Repl. Pad	Repl. Strap
Universal	EK087000	EK087400	EK087300

HCFA Code: L1832

Optional Shoe Insert

Shoe Size - Men	Shoe Size - Women	Shoe Insert Size	Left	Right
Up to 4 (23 cm)	Up to 5 (24 cm)	Small	SM007101	SM007201
5-9 (24-27 cm)	6-10 (24-28 cm)	Medium	SM007105	SM007205
10-14 (28-32 cm)	11-15 (29-32 cm)	Large	SM007109	SM007209

Post-Op Knee Bracing

Townsend ROM-ES



Flexion stop settings from 110 to 0 degrees in 10 degree increments. Patient friendly drop lock mechanism for fixating the brace at full extension. Locking hinge caps prevent the patient from re-adjusting stop settings. highly adjustable velcro closures
 A User-Friendly Design. Flexion and Extension Settings removable ends for shortening the brace are two unique features of Townsend post-operative ROM braces (knee brace models, only).

Product Features and Benefits

ROM-ES: 4 straps; full shell pads or cooler EZ pads; length adjustment from 12-25 inches; Drop locks fixate hinges at full extension; multiple extension and flexion stops (in 10 degree increments); Velcro tab strap closures or buckle closures (with full shell pads only); Synergistic Suspension Strap.

Size and Order Information

PART #	PAD TYPE	STRAP	LENGTH
TRESFPBS	FULL	BUCKLE	STANDARD
TRESFPVS	FULL	VELCRO	STANDARD
TRESEZVS	EZ	VELCRO	STANDARD
TRESFPBM	FULL	BUCKLE	MAXIMUM
TRESEZBM	EZ	BUCKLE	MAXIMUM
TRESFPVM	FULL	VELCRO	MAXIMUM